


Simple Roll-down Technique for Placement of the OB/Mobius®



1. Collapse the lower blue ring.



2. Insert it into the abdomen (first towards the head, then towards the feet).



Gently pull the yellow ring to ensure insertion and perform a quick digital exam to check for trapped tissue.



4. Position hands at 11 and 1 o'clock, using an assistant as the 12 o'clock position.



5. Roll the ring outwards and down as if scooping. If desired, an assistant can hold the rolled ring while physician slides fingers down the ring and completes the roll.



6. Repeat 2-3 times until desired exposure is achieved.



7. Perform procedure.



8. To remove, reach inside the abdomen.



9. Grasp the lower blue ring and pull outwards.



(203) 601-5200 | (800) 243-2974 | coopersurgical.com | 95 Corporate Dr. | Trumbull, CT 06611 Also on YouTube.com/coopersurgical