RUMI® II Handle

Reusable Uterine Manipulator Handle

UMH650

Directions for Use (English)

**NON-STERILE
STERILIZE BEFORE USE**

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

CE0086

DEVICE DESCRIPTION


The RUMI® II Uterine Manipulator Handle is a reusable device used with a sterile, single-use RUMI Tip.

The RUMI II Uterine Manipulator Handle (Handle) along with the RUMI Uterine Manipulator Tip (Tip) is used for positioning the uterus during laparoscopic surgery and in dye delivery where chromopertubation is necessary.

The Grip, Trigger, Arm and Snap Drum are constructed of medical grade autoclavable materials. The Snap Drum serves for attachment of the Tip.

Turning the Grip while depressing the Trigger permits the movement of the Snap Drum, thereby moving the position of the Tip. When the Trigger is released, it locks into the desired position.

Not made with natural rubber latex.

WARNINGS

- The uterus must be sounded for both depth and direction prior to application of the device.
- **DO NOT** use the Handle as a uterine sound.
- **NEVER** attempt uterine manipulation without a clear view of the uterus.
- As with all uterine manipulating devices, a careful clinical evaluation should be performed prior to use.
- Certain clinical conditions may present a uterus which is more prone to perforation and bleeding.
- Dye injection should be performed **SLOWLY**. Because of the efficient air/liquid seal created at the internal cervical os by the balloon, rapid injection of fluids (dye) may create intrauterine pressures, which could cause uterine damage and/or result in fallopian tube spasms.
- The Handle has only been tested for use with patients who have been anesthetized. Due to the need to dilate the cervix, the Handle is not recommended for use in non-anesthetized patients.

THE TIP IS SINGLE USE ONLY • NEVER REUSE A MANIPULATOR TIP

- Contents supplied sterile. Do not use if sterile barrier is damaged.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

INTENDED USE / INDICATIONS FOR USE

The Handle is designed for use in operative endoscopy (laparoscopy), where a uterus is present and where positioning of the uterus, fallopian tubes and ovaries or vagina is desirable. These types of surgeries include laparoscopic tubal ligation, diagnostic laparoscopy and/or operative laparoscopy. The Tips also provide for dye delivery in those procedures requiring chromopertubation.

CONTRAINDICATIONS

The Handle should not be used in patients who are pregnant or who are suspected of being pregnant, planned gamete intrafallopian transfer procedures, in patients who have an IUD in place, in patients with suspected pelvic infections and in cases where the surgeon deems it inadvisable or finds it difficult to insert the silicone tip into the cervix or uterus.

PRECAUTIONS

- Sterilize before use.
- Dilate the cervix to Hegar/Hank 8 (French 24) to ease patient insertion.
- Refer to the RUMI® Uterine Manipulator Tip Directions for Use for further precautions.
- Inspect Handle prior to use for proper operation.

ADVERSE REACTIONS


The following adverse reactions have been suspected or reported to be associated with all uterine manipulators. The order of listing does not indicate frequency or severity: cramping, infection, uterine and fallopian tube spasm with associative temporary physiological blockage of patient fallopian tubes and uterine perforation.

DIRECTIONS FOR USE

Note: Refer to the RUMI Uterine Manipulator Tips Directions for Use to ensure proper selection of Tip size and Tip attachment / detachment instructions.

1. Select a Tip which is less than or equal to the sounded depth of the uterus and attach it to the Handle.
2. Ensure the Tip Catheters are secure in the Catheter Channels of the Handle (see Figure 1).

Figure 1



UTERINE POSITIONING

The RUMI Handle is always intended to be used with the Grip in the upward position (see Figure 1).

TO ANTEVERT or insert into an anteverted uterus (see Figure 2): depress the Trigger and turn the Grip clockwise lifting the uterus to achieve the desired degree of anteversion, up to 90°. Release the Trigger to lock the desired position.

TO RETROVERT or insert into a retroverted uterus (see Figure 3): depress the Trigger and turn the Grip counterclockwise. A counterclockwise rotation of the Grip will result in the uterus being lowered to a retroverted position down to 50°. Release the Trigger to lock the desired position.

RECOMMENDED CLEANING INSTRUCTIONS

Before cleaning the Handle, remove the Tip as per the RUMI® Uterine Manipulator Tips Directions for Use.

1. Prepare the enzymatic detergent according to the manufacturer's recommendations and submerge the Handle in the detergent. Soak for ten (10) minutes.
2. Scrub the submerged Handle with a soft sponge and agitate. Clean lumens and crevices, using a pipe cleaner/brush or flush with a syringe when appropriate. Actuate any moving parts to loosen trapped soil.
3. Rinse in warm 100 °F - 120 °F (38 °C - 49 °C) tap water for one (1) minute. Thoroughly flush all lumens and other difficult to reach areas. Actuate while rinsing.
4. Ultrasonically clean the Handle for ten (10) minutes in a neutral pH detergent (Neutrad® or acceptable alternative). Prepare the detergent according to the manufacturer's recommendations.
5. Rinse the Handle with clean tap water, actuating moving parts while rinsing, for one (1) minute; repeat rinse twice.
6. Dry the Handle thoroughly with a clean, lint free cloth.
7. Repeat this cleaning procedure if the Handle appears to be soiled after initial cleaning and until there is no sign of soil on the device.

RECOMMENDED STEAM AUTOCLAVE STERILIZATION PARAMETERS

- The instrument must be thoroughly cleaned of all foreign matter prior to sterilization.
- Follow the manufacturer's instructions for operation and loading of steam autoclave.
- There must be direct steam exposure to all surfaces of the instrument being sterilized.

Pre-vacuum Cycle:

Preconditioning Pulses: 4

Exposure Temperature: 270 °F - 273 °F (132 °C - 134 °C)

Exposure time: 4 minutes

(continued)

Gravity Displacement Cycles:

Exposure Temperature: 270 °F (132 °C)

Exposure Time: 15 minutes

Exposure Temperature: 250 °F (121 °C)

Exposure Time: 30 minutes

WARRANTY

The RUMI® II Uterine Manipulator Handle is warranted by CooperSurgical for a period of 120 days from shipment. This warranty applies only to the original purchaser and only against defects in workmanship or materials, which under normal use render the instrument inoperable. The original purchaser shall prepay the shipping costs for returning a CooperSurgical instrument for warranty service. CooperSurgical will, at its option and without charge, either repair or replace any instrument, which CooperSurgical determines to be defective in material or workmanship. Instrument damage caused by misuse or accident, shall void this warranty.

CooperSurgical disclaims any liability for special, incidental, consequential, punitive or exemplary damages arising out of the use of any CooperSurgical instrument. CooperSurgical's liability in all events is limited to, and shall not exceed, the purchase price paid. Except as expressly provided in this limited warranty section, CooperSurgical makes no representations or warranties, expressed or implied, as to the instruments, including warranties of merchantability and fitness for a particular purpose, all of which are hereby expressly disclaimed and excluded. No warranty or affirmation of fact, express or implied, other than as set forth in this limited warranty, is made or authorized by CooperSurgical.

EXPLANATION OF SYMBOLS

 Reorder Number	 U.S. Federal law restricts this device to sale by or on the order of a physician.
 Serial Number	 Product conforms to the Medical Device Directive 93/42/EEC
 Caution	 Authorized Representative in the European Community.
 Consult instructions for use	 Manufacturer
 Do Not Reuse (Tips Only)	
 Keep Dry	

Not made with natural rubber latex.

Neutrad® is a registered trademark of Decon Labs, Inc.

All other marks bearing the symbols TM and ® are trademarks and registered trademarks of CooperSurgical, Inc.

© 2014 CooperSurgical, Inc. All rights reserved.

CooperSurgical

95 Corporate Drive • Trumbull, CT 06611 USA

Phone: (800) 243-2974

Fax: (800) 262-0105

www.coopersurgical.com

International

Phone: (203) 601-9818

Fax: (203) 601-4747

CooperSurgical, Inc.
95 Corporate Drive
Trumbull, CT 06611

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Neutrad® é uma marca registada da Decon Labs, Inc.
© 2014 CooperSurgical, Inc. Todos os direitos reservados.